

PERACETIC ACID 4.5% RENACIDIC 5 LTR

RENACIDIC PERACETIC ACID 4.5%

Main Properties


Effective against bacteria, fungi, spores and viruses, including hepatitis B viruses, hepatitis C viruses, HIV viruses.

Dissolves carbonate residues, especially suitable for disinfection and cleaning after bicarbonate dialysis.

Effective against the biofilm and effective against endotoxins (Pyrogens)

Easily rinsed off (Without salt out risk)

Produces no noxious odors in the diluted form



Composition

- Hydrogen Peroxide 21%
- Glacial Acetic Acid 15%

Application

Disinfection of:

- Dialyzers
- Medical apparatus
- Water treatment devices



Distributor by: